INSTRUCTIONS FOR USE

Insignis-26GTM

26 Gauge Subcutaneous Needle Sets

- EN** Insignis-26GTM Subcutaneous Needle Set – Instructions for use
- CS** Sada pro subkutánní podání (Insignis – 26GTM) - Návod k použití
- DA** (Insignis – 26GTM) Sæt med subkutane kanyler - Brugsanvisning
- DE** (Insignis – 26GTM) Subkutane Nadelsets - Gebrauchsanleitung
- EL** (Insignis – 26GTM) Σετ βελονών υποδόριας χορήγησης – Οδηγίες χρήσης
- ES** Equipos de agujas subcutáneas (Insignis – 26GTM) – Instrucciones de uso
- FI** (Insignis – 26GTM) Subkutaaniset neulasarjat – Käyttöohjeet
- FR** Kits d'aiguilles sous-cutanées (Insignis – 26GTM) – Mode d'emploi
- HE** הוראות שימוש (Insignis – 26GTM) ערכות מחט תת עורית
- IT** Set di aghi sottocutanei (Insignis – 26GTM) – Istruzioni per l'uso
- NL** (Insignis – 26GTM) subcutane naaldsets – Gebruiksaanwijzing
- NO** (Insignis – 26GTM) subkutane nålesett – Bruksanvisning
- PL** (Insignis – 26GTM) Zestawy igłowe do podawania podskórnego – Instrukcja użytkowania
- PT** Conjuntos de Agulhas Subcutâneas (Insignis – 26GTM) – Instruções de utilização
- RO** (Insignis – 26GTM) Seturi de ace subcutanate – Instrucțiuni pentru utilizare
- SK** Súprava subkutánných ihl (Insignis – 26GTM) – Návod na použitie
- SV** (Insignis – 26GTM) Subkutana nålset – Bruksanvisning
- TR** (Insignis – 26GTM) Subkütan İğne Setleri – Kullanım talimatı
- UK** Набір голок для підшкірного введення Insignis-26GTM – Інструкції із застосування
- RU** Набор игл для подкожного введения Insignis-26GTM – инструкции по применению

INNOVATIVE HEALTH SCIENCES LLC


1108 Kings Highway, Suite #4
Chester, NY 10918 USA
+1-855-680-0630

www.innohealthsci.com

04-1027 Rev 05 2026-08-14

Insignis-26G™ Subcutaneous Needle Sets (box of 10)

Configuration	6mm	9mm	12mm	14mm
1 site	SC-NS-0106	SC-NS-0109	SC-NS-0112	SC-NS-0114
2 sites	SC-NS-0206	SC-NS-0209	SC-NS-0212	SC-NS-0214
3 sites	SC-NS-0306	SC-NS-0309	SC-NS-0312	SC-NS-0314
4 sites	SC-NS-0406	SC-NS-0409	SC-NS-0412	SC-NS-0414

Sterile dressing/adhesive is included in the administration set box.

	Caution	STERILE R	Sterilized using irradiation	Rx Only	Prescription only
	Consult instructions for use		Do not resterilize	MD	Medical device
	Do not use if package is damaged		Do not reuse		Sterile barrier system
	MR Unsafe		Temperature limit		Humidity limitation

EN Indications for use: The Insignis-26G™ Subcutaneous Needle Sets are used for the subcutaneous infusion to administer fluids from a container or a syringe to a patient's subcutaneous tissues.

Maximum duration of use: < 24 hours

**Caution**


- Use of the Insignis-26G™ Subcutaneous Needle Sets are intended only for the patient for whom the device is prescribed.
- U.S. federal law restricts this device to sale by or on the order of a physician.
- Before use, carefully inspect the administration set package. Do not use if the package has been opened or is damaged.
- The administration sets are single use only. Reusing a single-use device could cause the device to not perform as intended or may result in infection.
- Do not re-sterilize administration sets.
- It is recommended to store the administration sets in a dry place (humidity levels 15–90%) at room temperature (20–25°C; 68–77°F).
- Do not stretch or kink the administration sets/ancillary devices.
- Always follow the infusion instructions located on the package insert of the manufactured drug.
- The Insignis™ Syringe Infusion System has not been tested in a Magnetic Resonance Imaging (MRI) environment and therefore cannot be recommended for use in such an environment.
- If using the infusions set for the first time, do so in the presence of a healthcare provider.
- Since individuals inherently have varying amounts of subcutaneous fat on the body, be sure to select the appropriate needle length to fit your needs or as instructed by a physician.
- If the infusion site becomes inflamed, replace the set, and choose a new infusion site until the first site has healed.
- Be sure to remove the needle guard prior to insertion.
- Ensure that the needle is not bent prior to insertion.
- Never attempt to free a clogged tubing line while the set is inserted.
- Store infusion sets in cool, dry place. Do not leave infusion sets in direct sunlight or inside a vehicle.
- When using multi-leg needle set configurations, over-infusion may occur to remaining sites if one or more of the needles is blocked.
- Inaccurate medication delivery, infection, and/or site irritation may result from improper needle insertion/maintenance of infusion site.
- Always consult and comply with guidance provided by your healthcare providers, and drug manufacturers, regarding the location and the number of infusion sites.

Adverse Event Reporting


Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Contact information for both is shown on the last page of this document.

EN Starting the Infusion for Subcutaneous Administration

1. Wash hands thoroughly and if needed, wear disposable gloves.
2. Verify that you have selected the appropriate needle set for your application.
3. Remove end cap and connect to the infusion source.
4. Using aseptic technique, prep infusion site(s) as directed by a clinician.
5. Remove rubber band from butterfly and fold wings back to expose needle guard; remove needle guard.
6. Insert needle(s) dry into appropriate locations in subcutaneous tissue at 90° angle.
7. Lay butterfly wings flat on skin. Cover infusion site with needle dressing/adhesive tape. Start infusion as recommended by your provider.

**Ending the Infusion**

1. Pause before removing needles; allow site pressures to decrease naturally. Remove needle dressing. To minimize accidental risk of needle stick injury, use one hand to remove needle. Using the same hand, squeeze butterfly safety closure wings together to firmly encase needle.
2. Once wings are closed, dispose of needle set in sharps container.

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
Definition of Symbols

	Caution	REF	Catalog number	SN	Serial number
	Consult instructions for use	LOT	Batch code	MD	Medical device
	Use by YYYY-MM-DD	QTY	Quantity	UDI	Unique device identifier
	Date of manufacture	STERILE R	Sterilized using irradiation	Rx Only	Prescription only
	Temperature limit		Do not resterilize		Do not reuse
	MR Unsafe		Humidity limitation		Do not use if package is damaged
	Sterile barrier system	EU REP	European Authorized Representative		Manufacturer
CE	European Conformity	CH REP	Swiss Authorized Representative		Importer

OncoMedical AG
 Hintere Hauptgasse 9
 CH-4800 Zofingen,
 Switzerland

MedEnvoy Global B.V.
 Prinses Margrietplantsoen 33 –
 Suite 123 2595 AM The Hague
 The Netherlands

UKRP

MedEnvoy UK Limited
 167–169 Great Portland Street
 5th Floor
 London, W1W 5PF
 United Kingdom

INNOVATIVE HEALTH SCIENCES LLC
 1108 Kings Highway, Suite #4
 Chester, NY 10918 USA
 +1-855-680-0630
www.innohealthsci.com

Insignis-26G™ Subcutaneous Needle Sets is a trademark of Innovative Health Sciences LLC (IHS). IHS is ISO 13485 certified and abides by the quality policies as outlined in 21 CFR 820.